


Fit Your Needs,  
**GM85 Fit**-type




# FOCUS ON THE CORE OF "MOBILE" RADIOGRAPHY

In recent extraordinary situations, radiography is facing a new paradigm in its ways of working. As the procedure volume of mobile x-ray increases to deal with unprecedented patient volume, Samsung is focusing more on efficient mobile imaging based on its accumulated technology. Experience higher productivity and performance with GM85 Fit designed to fit your workflow, diagnosis and care.



- Higher efficiency
- Shorter service time
- Longer lifespan



<p><b>Rapid Charging<sup>1</sup></b></p> <p>Within <b>3~4</b> hours</p>	<p><b>All Day Operation<sup>2</sup></b></p> <p><b>500</b> exposures</p> <p><b>50</b> km (31 miles)</p> <p><b>10</b> hours</p>	<p><b>Max. Exposure<sup>3</sup></b></p> <p><b>2,300</b> exposures</p>	<p><b>Sleep Mode Stand-by<sup>4</sup></b></p> <p><b>31</b> hours</p>
---	---	---	--

1) 200 ~ 240 VAC < 3 hours, 100 ~ 127 VAC < 4 hours  
 2) Test condition : Chest AP / 80 kVp / 160 mA / 5 msec / 60 sec intervals, with driving at 5.6 km/h  
 3) Test condition : Chest AP / 80 kVp / 160 mA / 5 msec / 10 sec intervals, without driving  
 4) Stand-by 18 hours for LCD on, 31 hours for sleep mode, and 75 hours for power off

# Fit Your Workflow

WHEN YOU GO TO START AN EXAM,




At the moment you start working, GM85 Fit increases your work efficiency with the fast speed of booting up the system. **Quick Driving** allows radiologists to move GM85 fit right away without waiting for the system to be fully ready. Users can easily log in to the system through the **embedded RFID badge tag\***. It is a built-in type in the equipment, so it reduces external interference.

**Additional Patient Information\*** is useful when checking the list of patients requiring x-ray imaging. It intuitively shows infection and fall risk information and orders comments that support correct action for each patient. This helps keep patients and staff safe by responding to possible risks and taking preemptive measures. The **lightweight** GM85 Fit allows easy maneuvering both in and out of elevators without worrying about the weight limit.

Embedded RFID badge tag\*




Additional Patient Information\*



# Fit Your Workflow

## WHEN YOU REACH A PATIENT,




## WHEN YOU GET X-RAY IMAGES,



When checking the patient's identity after arriving at the room, a **wireless barcode scanner\*** allows you to freely approach the patient and scan information without being restricted by obstacles that may be near you and the patient. With a **wireless hand switch** that supports up to 10 m (33 ft) range, users can take images from far away to reduce radiation exposure.


### Prior Exam Review\*



Since the role of mobile DR is usually to use follow-up exams, it is important to perform x-ray exams under the same conditions. **Prior Exam Review\*** displays previous images and exposure parameters of the patient being examined. With a quick comparison, users can improve image consistency and reduce retakes.

**SimGrid™\*** application streamlines the workflow by guaranteeing image quality without the use of a conventional grid. This allows the omission of grid installation and removal steps from the conventional workflow leading to 28% reduction in total exam time.

### SimGrid™\*



In portable exams, radiologic technologists rotate images after the acquisition as the image orientation depends on the exam environment. **Image Auto rotation\*1** detects the rotated angle of the image and automatically rotates in the correct direction based on AI algorithm. (97% average accuracy) The **auto-run option** of advanced applications creates the companion image without additional settings or x-ray exposure.

In the ER/OR or Trauma environments, multiple staff may often need to quickly acquire and check images in real-time. **Mirror View\*** provides secured screen sharing using Wi-Fi CERTIFIED Miracast™. This allows staff to check images together on a separate screen, reducing the time to first aid and the risk of contamination. Through **QuickLink**, you can check whether the acquired images uploaded well to the server by accessing RIS/HIS directly without a separate PC.

When preparing for the next exam, **EasyClean**, which locks the screen for a while, helps you clean equipment quickly and control infection effectively.

### Water resistant detector for easy cleaning




### Detector rest for easy bagging



1) [Chest/Abdomen/Pelvis/Infantogram] 0°, 90°, 180°, 270° Rotation  
[Hand/Foot] 1° Unit Rotation

# Fit Your Diagnosis



Up to **50%** Dose reduction  
Same image quality

Samsung's S-Vue™ image processing engine provides excellent, high resolution images for diagnostic accuracy and confidence. The S-Vue™'s adaptive filtering and advanced denoising technology offer the natural presentation of bone and soft tissues. Also, dose level can be reduced up to 50% for adult chest and 47.5% for adult abdomen without compromise in image quality.

Note : The claim concerning Samsung DR is based on limited phantom and clinical study results. Only routine PA chest radiography and abdominal radiography for average adults and pediatric abdominal, chest, skull radiography were studied, excluding pediatric patients under 1 month old. (FDA cleared - K172229, K182183) In practice, the values of dose reduction may vary accordingly. These clinical images calculates the dose reduction rate from its own standard dose at the clinical site, unlike our FDA claim which compares dose between new IPE and old IPE. The clinical site is responsible for determining whether the particular radiographic imaging needs are not impacted by such x-ray dose reduction.


**ACC** Standard Detector  
ACCESS · ACCURACY · EFFICIENCY

IP54




High quality  
Small pixel pitch

To offer optimized image quality, Samsung introduce AccE Standard Detector that broaden customers' choice of Csi detectors. This latest lineup will meet various clinical needs for sizes. Enhanced load allowance along with dust and water resistance allows the detector to be actively implemented in versatile environments such as ER and OR.



1) Based on 4cm radius on the center

Case 1. Chest AP (X-Large adult)



The image was taken with GM85.

Case 2. Abdomen



The image was taken with GC85A.

Case 3. L-spine



The image was taken with GC85A.

Case 4. Pelvis




The image was taken with GC85A.

## Consistent Image Quality

The upgraded S-Vue™ displays areas with varying bone density clearly regardless of the patient's size or body position through advanced region-specialized processing. With this engine, halo artifacts surrounding metallic objects are mitigated and visualized properly.

# Fit Your Diagnosis




## On-Device CAD\*

Samsung's new On-Device CAD solution with VUNO Med®-Chest X-ray™ assists you to diagnose more efficiently and accurately by automatically detecting the 5 most common thoracic findings in emergency care.


- 5** Nodule/Mass · Consolidation · Interstitial Opacity
- Pleural Effusion · Pneumothorax

Case. Chest AP

Images were taken with GM85.



Without On-Device CAD



With On-Device CAD




## SimGrid™\*

With just a click, SimGrid™ allows you to provide better patient care with higher satisfaction and reduced retake rates without the use of a portable grid. It improves image contrast by reducing scatter radiation effects and creates better image quality. The 3-step intensity control (Low/Medium/High) enables customized image processing.

Case. Chest AP

Images were taken with GM85.



Without Grid



With SimGrid™




## S-Enhance\*

To support your diagnosis, S-Enhance improves the clarity of foreign bodies (e.g. tube, line and/or needle) in images of chest, abdomen, and L-spine. With a single on-screen click, the companion image is created without additional settings or x-ray exposure, streamlining the workflow.

Case. Chest AP

Images were taken with GM85.



Without S-Enhance



With S-Enhance



## Bone Suppression\*


Without additional setting or exposure, Bone Suppression Imaging improves the clarity of soft tissues by suppressing the appearance of bones in chest images, which improves your ability to detect nodules. You can easily create the companion image with just a click on the screen.

Case. Chest PA

Images were taken with GC85A.



Without Bone Suppression



With Bone Suppression

# Fit Your Care




Up to  
**45%**

## Dose reduction For very young patients


Underaged patients are more radiosensitive than adults. To alleviate these concerns, the new S-Vue™ engine helps achieve the optimal dose level for children during pediatric x-ray scans. The dose level can be reduced up to 45% for pediatric abdomen, 15.5% for pediatric chest, and 27% for pediatric skull exams. This is especially significant as abdomen protocols may include genital regions.

### Case. Pediatric Chest AP

Images were taken with GM85.



**15%**  
DOSE  
REDUCTION




( 54 kVp / 1.42 mAs / 0.06 dGy\*cm<sup>2</sup> / 0.1 mmCu Filter)

( 54 kVp / 1.21 mAs / 0.05 dGy\*cm<sup>2</sup> / 0.1 mmCu Filter)

## 6stages Pediatric Exposure Management

Optimized **6-stage weight dependent imaging** enables pediatric patients to avoid unnecessary x-ray exposure using precise dose management, resulting in superior image quality. Users who need detailed control can apply 3 types of **additional filters**. (0.1/0.2/0.3 mm Cu)



## Less than 43dB when the system is on stand-by

GM85 Fit allows low noise operations that do not disturb other patients. With **night mode** on, the screen color, temperature, and sound volume are accordingly adjusted. It will be a necessary function for taking x-ray exams in a quiet environment such as NICU. Its **quick exposure** feature is useful for pediatric patients who are not easy to maintain stationary posture.



American Academy of Pediatrics (AAP) recommendation in NICU ≤ 45 dB



## Efforts to reduce fatigue of radiographers

The **lightweight** GM85 Fit provides an effortless driving experience for users. Not only the equipment but also the detector that must be held directly with the arm should be light. The user-centric design of AccE Standard Detector supports patient positioning and alleviates daily burdens.

**339 kg**  
(747 lb)  
Light weight

**CENTER ENGRAVING**  
to help position the patient



**REAR GRIP**  
to support transportation

**SIDE CHAMFER**  
to ease your lifting



## Cybersecurity & Value Care Service

Samsung provides a solution to support our customers by offering the tools to protect against cyber threats that may compromise invaluable patient data and ultimately degrade the quality of care. In addition, **remote software update** function allows the system to keep up to date. Experience Samsung's **value care service** that improves the system's uptime and total cost of ownership with 26% faster repair service.<sup>1</sup>



<sup>1</sup>) Based on the service time evaluation for major parts repair (compared to the conventional GM85)

About Samsung Electronics Co., Ltd.

Samsung Electronics Co., Ltd. inspires the world and shapes the future with transformative ideas and technologies.

The company is redefining the worlds of TVs, smartphones, wearable devices, tablets, cameras, digital appliances, medical equipment, network systems, and semiconductor and LED solutions. For the latest news, please visit the Samsung Newsroom at [news.samsung.com](https://news.samsung.com).

### **GM85 Fit Catalog v1.2-230217-CE**

GM85 Fit is a new configuration of GM85. (Fit-type)



Scan code or visit  
[www.samsunghealthcare.com](https://www.samsunghealthcare.com)  
to learn more

**Represented by:**  
**Commonwealth X-Ray, Inc.**  
**[www.commonwealthxray.net](https://www.commonwealthxray.net)**  
**[cxr@cxrinc.net](mailto:cxr@cxrinc.net)**  
**859-885-4854**

Copyright © 2022 Samsung Electronics Co. Ltd. All rights reserved. Samsung is a registered trademark of Samsung Electronics Co. Ltd. Specifications and designs are subject to change without notice. Non-metric weights and measurements are approximate. All data were deemed correct at time of creation. Samsung is not liable for errors or omissions. All brand, product, service names and logos are trademarks and/or registered trademarks of their respective owners and are hereby recognized and acknowledged. Samsung Electronics Co., Ltd.

129, Samsung-ro, Yeongtong-gu, Suwon-si, Gyeonggi-do, 16677, Republic of Korea

**SAMSUNG**